

S131 PIONEER 3 3-WHEEL
MIDI ELECTRIC SCOOTER



Order Form (Effective May 1, 2017)

For ease of ordering contact Inside Sales toll free at 1-800-963-7487, Fax: 239-574-2661

Order Taken by: _____

Company Name: _____ Account # _____

Phone: _____ Date: _____ P.O # _____

Ship to Name/Address: _____

Special Note:

DME suppliers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.

SCOOTER TYPE	PRICE	HCPCS
<input type="checkbox"/> S131 Pioneer 3 (RED).....	\$3198.00	K0807
<input type="checkbox"/> S131 Pioneer 3 (BLUE).....	\$3198.00	K0807










Battery

<input type="checkbox"/> U1 Battery (each).....	\$158.00
---	----------

OPTIONS & ACCESSORIES

<input type="checkbox"/> Oxygen Tank Holder	\$138.00
<input type="checkbox"/> Cane & Crutch Holder.....	\$98.00
<input type="checkbox"/> Walker Holder.....	\$118.00
<input type="checkbox"/> Foldable Rear Basket (Med.).....	\$72.00
<input type="checkbox"/> Foldable Rear Basket (Lrg.).....	\$108.00
<input type="checkbox"/> Cover.....	\$180.00
<input type="checkbox"/> Canopy.....	\$298.00

DIMENSIONS

 length	46	 seat depth	20
 width	23	 back height	20
 height	45 1/2	 seat height from deck	16-19
 wheel base	32	 seat height from ground	25-28
 seat width	20		

PRODUCT SPECIFICATIONS

Max Speed:	5 mph
Range up to:	18 miles
Turning Radius:	31"
Ground Clearance:	4"
Motor:	DC 24V, 250W
Controller:	Dynamic Rhino 70A
Battery:	12V/ U1 (34AH)
Charger:	5A off-board
Gradient:	8°
Tires	
Front Wheel	10" pneumatic tire
Rear Wheels	10" pneumatic tires
Brake:	intelligent, regenerative electromagnetic brakes
Product Weight	
Base	86 lbs
Battery (each)	26 lbs
Seat	40 lbs
Total Weight with Battery	178 lbs
Weight Capacity	450 lbs
Warranty	
Frame	Limited Lifetime
Electronics	1 year
Motor/Brakes	1 year
Batteries	1 year

Features:

- High-torque drive train allows for increased power and maneuverability
- Adjustable tiller to fit various arm lengths and collapses for easy storage
- Dashboard control panel for ease of use
- Flip-up armrests with height, width & angle adjustment
- Steel bumper protects shroud from scratches and serves as a lift handle
- One Piece frame design
- Limited lifetime warranty on frame
- Sliding captains seat



Merits Health Products, Inc.,
730 NE 19th Place,
Cape Coral, FL 33909
www.meritshealth.com

Quantity for above selection