

**S235 PIONEER 1 3-WHEEL
MINI ELECTRIC SCOOTER**



Order Form (Effective May 1, 2017)

For ease of ordering contact Inside Sales toll free at 1-800-963-7487, Fax: 239-574-2661

Order Taken by: _____

Company Name: _____ Account # _____

Phone: _____ Date: _____ P.O # _____

Ship to Name/Address: _____

Special Note:

DME suppliers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.

SCOOTER TYPE	PRICE	HCPCS
<input type="checkbox"/> S235 Pioneer 1 (RED).....	\$2198.00	K0806
<input type="checkbox"/> S235 Pioneer 1 (BLUE).....	\$2198.00	K0806










PRODUCT SPECIFICATIONS

Max Speed:	4 mph
Range up to:	15 miles
Turning Radius:	27"
Ground Clearance:	3 ½"
Motor:	210W
Controller:	Dynamic Rhino 70A
Battery:	12V/22AH x 2
Charger:	3A off-board
Gradient:	10°
Tires	
Front Wheels	9" foam filled
Rear Wheels	9" foam filled
Brakes:	intelligent, regenerative electromagnetic brakes
Product Weight	
Base	73 lbs
Battery (each)	14 ½ lbs
Seat	27 ½ lbs
Total Weight with Battery	130 lbs
Weight Capacity	300 lbs
Warranty	
Frame	1 year
Electronics	1 year
Motor/Brakes	1 year
Batteries	1 year

OPTIONS & ACCESSORIES

<input type="checkbox"/> Oxygen Tank Holder	\$138.00
<input type="checkbox"/> Cane & Crutch Holder.....	\$98.00
<input type="checkbox"/> Walker Holder.....	\$118.00
<input type="checkbox"/> Foldable Rear Basket (Med.).....	\$72.00
<input type="checkbox"/> Foldable Rear Basket (Lrg.).....	\$108.00
<input type="checkbox"/> Cover.....	\$180.00
<input type="checkbox"/> Canopy.....	\$298.00

DIMENSIONS

 length	44	 seat depth	16
 width	23	 back height	15
 height	34.5	 seat height from deck	16-20
 wheel base	30	 seat height from ground	19.5-23.5
 seat width	17		



Merits Health Products, Inc.,
730 NE 19th Place,
Cape Coral, FL 33909
www.meritshealth.com

Quantity for above selection

Features:

- Compact design and easy disassembly
- Comfortable adjustable seat and tiller to fit various individual users.
- Tight turning radius provides excellent maneuverability
- Ideal for both indoor and outdoor use
- Height and width adjustable armrests
- Adjustable tiller